
Juno Targeted DNA Sequencing Library Preparation System

Affordable scalability for mid-to-highthroughput next-generation sequencing





Scale your targeted DNA library preparation using Juno and multiple sample throughput options, up to 192 samples per processing run and thousands of genomic loci targeted per sample. Design your custom NGS panel with ease, either online or with our consultative, expert assay design services. You may add amplicon assays over time so your panel content will always be up to date.

Highlights

Easily scale library preparation from dozens to thousands of samples per week.

Control costs by reducing reaction volumes and total hands-on time per sample.

Design with ease your custom NGS panels from a few to thousands of genomic variants.

Targeted DNA Sequencing Library Preparation	LP 48.48 IFC	LP 192.24 IFC
Number of samples per run	Up to 48	Up to 192
Number of amplicon assays per sample	Up to 4,800	Up to 2,400
Amplicon length	150–500 base pairs	150–500 base pairs
Workflow time from DNA to NGS-ready libraries (1 run)	8 hours for 48 NGS-ready libraries, including 4 hours of walkaway automation	9.5 hours for 192 NGS-ready libraries, including 5.5 hours of walkaway automation
Sample throughput in 24 hours	Up to 96 sample libraries	Up to 384 sample libraries
Number of samples multiplexed per NGS run	Up to 1,536 using Fluidigm barcodes	Up to 1,536 using Fluidigm barcodes

Juno Library Preparation Workflow

1

2

3

4

DESIGN PANEL

Enter genes or genomic coordinates online using the D3™ assay design portal.

TRANSFER

Transfer DNA samples, assay primers and other reagents to the integrated fluidic circuit (IFC). Samples and reagents combine controlled reactions via the IFC.

CONSTRUCT LIBRARIES

Automate workflow for targeted amplification and sample barcoding within the IFC. Adapter addition and cleanup are completed in a single tube.

SEQUENCE

Perform NGS using an Illumina® system. With these libraries you have the capability of multiplexing up to 1,536 samples per sequencing run.

Maximize high-throughput NGS

Scale your NGS library production from dozens to thousands of libraries per week. Choose from multiple sample formats to meet your lab's needs for throughput and turnaround time.

Realize operational efficiency

Lower your reagent and primer costs by reducing individual reactions to nanoscale volumes using microfluidic technology. Achieve additional labor savings by reducing total hands-on time per sample in mid-to-high-throughput settings.

Simplify NGS panel design

Our D3 assay design delivers expert design support and optimized PCR primers in a format tailored to your needs.

Ordering information

Product Name	Part Number	IFC Type and Quantity
Targeted DNA Seq Library Assays	ASY-MPX ASY-MPX-P	Various
Targeted DNA Seq Library Adapter Set	101-2412	20 IFCs
Targeted DNA Seq Library Barcode Plates	101-0744	192.24–40 IFCs 48.48–160 IFCs
Advanta™ NGS Library Prep Reagent Kit	101-7663 101-7666 101-7667 101-7669	LP 48.48, 2 IFCs LP 48.48, 10 IFCs LP 192.24, 2 IFCs LP 192.24, 10 IFCs
Juno Accessories Kit (IFCs with accessories)	101-2334 101-2333 101-2349 101-2348	192.24–5 IFCs 192.24–10 IFCs 48.48–5 IFCs 48.48–10 IFC

Learn more at

fluidigm.com/juno

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Patent and license information:** fluidigm.com/legalnotices. Fluidigm, the Fluidigm logo, Advanta, D3 and Juno are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2018 Fluidigm Corporation. All rights reserved. 07/2018

